

Batch Mode

Codonics SLS 500i Point of Care Station (PCS) is perfectly suited for pharmacy or the hectic environments of the OR, ICU and PACU where duplicate labels are often required.

In the hospital pharmacy, it is commonplace to prepare medications in batches where multiple copies of the same label are required. The OR, ICU and PACU often have the need for multiple identical labels for IV applications, where the IV bag, the line, the catheter, and the pump must all be labeled alike. Multiple copies of the same label is a requirement when multiple syringes from a single multi-dose vial are prepared in advance or if multiple syringes of the same medication are required for a single case.

Based on the hospital's on-board, site-specific formulary, the system electronically verifies each and every vial and/or ampoule in hand as a first line of defense to prevent medication errors, such as with dangerous look-alike drugs and hard to read vials and ampoules.

With a click from the touch screen to select "batch" printing, a user simply enters the number of identical labels to be printed and, if desired, the lot number. When the user is ready, they quickly scan the container barcode of the drug in hand on SLS PCS and the system visually and audibly identifies the drug and concentration. SLS PCS prints easy-to-read, full color labels on demand that contain all of the elements required by the Joint Commission, and labels can be immediately applied to syringes. The user can select to have the labels presented individually or in a strip. Labels can also accommodate two sets of initials as needed, when one person prepared the label but not the syringe or when a second check for safety is required.

SLS PCS showing
Batch Mode

 Fentanyl
10 mcg/mL



KEY BENEFITS

Additional Advantages:

- Electronically identifies parenteral vial to prevent drug selection errors
- Fast, easy touch screen entry of the desired number of identical labels to be printed (0 to 999)
- Optional on-screen entry of lot number for each batch of labels to be printed
- Touch screen selection of individually cut or strip cut labels
- Touch screen display shows printing status (such as "Printing 10 of 250")
- Configurable per unit



Optional lot number



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
info@codonics.com
www.codonics.com

Call Codonics today at **800.444.1198** or visit
www.codonics.com for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: www.codonics.com/ip/patents Printed in the U.S.A. Copyright © 2016 Codonics Inc. 9/2016



Sponsoring Donor



Corporate Member