

# Codonics® SLS 500i® Operating Software

## Release Notes

### Version 2.1.3 SLS 500i Release Notes

#### Summary

Version 2.1.3 software is a software release for the Codonics Safe Label System® 500i (SLS) Point of Care Station (PCS) and replaces 2.1.1 software. SLS PCS does not support software downgrades to previous software versions.

**NOTE:** The 2.1.0 Administration Tool (AT) software must be used with the 2.1.3 SLS PCS software.

The 2.1.0 AT works with all 2.1.x SLS PCSs simultaneously.

As of this release, SLS software version 2.1.3 version is considered the most current, stable, safe and effective version for the **updated** SLS PCS. SLS software version 2.1.2 is considered the most current, stable, safe and effective version for **previous** SLS PCS. Codonics recommends and encourages updating to these versions.

**NOTE:** The serial number for the **updated** SLS PCS with the new printer begins with the number 142C.

The **previous** SLS PCS serial number began with either the number 140C or 141C.


**NOTE:** 2.1.3 SLS SW will not install on **previous** versions of SLS PCS (i.e., those starting with 140C and 141C).

This document provides information about Product Corrections (i.e., Corrective Updates). The SLS PCS User's Manual version 1.3.0 and User's Manual Addendum 1.8.0 provide further information about SLS PCS.


**WARNING:** The use of label colors is intended only as an aid in the identification of drug groups and does not absolve the user from the duty of reading the label and correctly identifying the drug prior to use.

#### Product Corrections

- ◆ **Multi-package updates can now be installed on SLS PCS.** SLS Product Bulletin 406 identified an issue where if a multi-package update was installed on SLS PCS, the wrong drug could be displayed on SLS. This issue no longer occurs. Refer to SLS Release Notes 901-245-014 and SLS AT File Manager Technical Brief 901-276-001 for details on multi-package updates.



- ◆ **Previously using large labels (80mm x 30 mm) with the updated SLS PCS, there could be smearing on the label.** During the last third of a roll, the labels can curl to the point where the last part printed on the label could smear. This will no longer occur on SLS PCSs that have either a Mod 7 or a Configuration 003 and higher without Mod 7. See serial number label for details.



- ◆ **Previously installing a 2.1.0 APP Package on a higher version of SLS SW (e.g., 2.1.1) using a USB drive would hang-up during installation.** Now, the SLS will not allow a 2.1.0 APP package to be installed using a USB drive on a SLS PCS that has a higher version (i.e. 2.1.2 or 2.1.3) of SLS SW. It will display the 2.1.0 APP package as INVALID.

## **Known Common Issues**

This section details common issues with updated SLS PCSs that are likely to be experienced by most users. These issues are specific to the current release; previous issues can be found in previous Release Notes. None of the issues included in the document are considered to render the product unsafe or ineffective for clinical use.

- ◆ **Pressing the Load button does not cut labels on the updated SLS PCS.** When pressing the Load button, the unit will feed media forward but does not cut it when the door is closed. Once a subsequent syringe label is printed, the unit will eject the printed syringe label attached to the blank label from pressing the Load button. When the Load button is pressed on previous SLS PCSs with serial numbers 140C and 141C, the unit will feed media forward. When the door is closed, it will cut a label.

## **Known Uncommon Issues**

This section details uncommon issues with updated SLS PCSs that are unlikely to be experienced by most users. These issues are specific to the current release; previous issues can be found in previous Release Notes. None of the issues included in the document are considered to render the product unsafe or ineffective for clinical use.

- ◆ **SLS PCS may not immediately connect to 5 GHz WiFi.**  
The SLS will report Offline in the AT-Device Manager but after 5 minutes it will be connected and report as Ready.

## **Technical Support**

If problems occur during software installation or operation, contact Codonics Technical Support at any time.

Phone: +1.440.243.1198

Email: [support@codonics.com](mailto:support@codonics.com)

Website: [www.codonics.com](http://www.codonics.com)

---

*Get it all with just one call*  
**1.800.444.1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2010-2016 Codonics, Inc. Printed in the U.S.A. Part No. 901-243-023.02



17991 Englewood Drive  
Cleveland, OH 44130 USA  
+1 (440) 243-1198  
+1 (440) 243-1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)

Codonics Trading Co, Ltd.  
317 Xianxia Rd. Building B  
Unit 1412  
Changning Dist., Shanghai  
P.R. China, 200051  
86-21-62787701  
86-21-62787719 Fax

Codonics Limited KK  
AQUACITY 9F,  
4-16-23, Shibaura  
Minato-ku, Tokyo, 108-0023  
JAPAN  
81-3-5730-2297  
81-3-5730-2295 Fax