

# Codonics® Safe Label System® Administration Tool Software

**Release Notes**

## **Version 2.6.1 Safe Label System Administration Tool Release Notes**

### **Summary**

Version 2.6.1 is a software release for the Codonics Safe Label System® (SLS) Administration Tool (AT). The software is compatible with SLS Point of Care Stations (PCS) installed with 2.6.1 software. As of this release, this version of software is considered the most current, stable, safe and effective version for SLS PCS.

This document provides information about Product Corrections (i.e., Corrective Updates) in 2.6.1 SLS AT software.

**NOTE:** The AT Server application provides all the same capabilities as the existing AT Stand-alone application. As of 2.5.0, the AT Stand-alone application is no longer supported.

**WARNING:** If you're using WIN 10 as your computer's operating system to run the AT, pay careful attention to uninstalling previous versions (prior to 2.5.0) of the AT. WIN 10 can inadvertently uninstall the current version (i.e., 2.5.0) of the AT. Contact Codonics Tech Support with questions about the correct process.

The SLS AT User's Manual version 1.3.0 and User's Manual Addendum 1.8.0 provides further information about the AT.

### **Recommended - Requirements for AT Server:**

**Processor:** Xeon W-2135 @ 3.70GHz, Core i7-7700 @ 3.60GHz or faster.

**RAM:** 8 GB for application.

**Application Disk Space:** 30 GB total at install.

### **Minimum - Requirements for AT Server:**

**Processor:** Xeon X5660 @ 2.80 GHz, Core i7-2600 @ 3.40GHz or faster.

**RAM:** 6GB for application.

**Application Disk Space:** 10 GB total at install

### **Additional Information for AT Server:**

**Computer OS:** Windows® 10, Windows Server 2012R2, 2016 and 2019. 64-bit required.

**Concurrent access:** Not supported.

### **Notes:**

- ◆ For best performance hosting the SLS AT data directory, use a locally attached SSD or high speed RAID.
- ◆ When running the SLS AT on a VM, allocate at least two CPU cores to the VM.
- ◆ When running more than one instance of the SLS AT on the same hosting system, increase the system RAM by at least 3 GB for each additional instance.
- ◆ Both configurations shown are sufficient for running the SLS Email Notifier (EN) on the same system hosting the SLS AT.

### **Recommended requirements for AT Client:**

**Browser support:** Internet Explorer IE 11 or newer, Edge, Chrome

**Screen Resolution:** Minimum 1440 x 900; use default browser fonts

**Computer OS:** Windows 10 or newer

**Applications:** Excel 2003 or greater to open reports

### **Product Corrections**

- ◆ **Certificates with multiple Subject Alternative Names (i.e., SAN) previously would prevent AT servers from starting.** In 2.6.0 SW, the AT would be able to import certificates using the AT server manager, but the AT server manager would fail to start. This is now fixed.

### **Known Common Issues**

This section details common issues with 2.6.1 AT software that are likely to be experienced by most users. These issues are specific to the current release; previous issues can be found in previous Release Notes. None of the issues included in the document are considered to render the product unsafe or ineffective for clinical use.

- ◆ **WARNING: The formulary report does not sort duplicates.** Duplicate Container IDs in the formulary are reported but are displayed in a way that may not put the duplicates together. You can sort the Excel spreadsheet on Container ID.

- ◆ **WARNING: Replacing drug names with non-Latin characters (e.g., Kanji, Hebrew) causes duplicate drug entries in the MDD.** When the MDD screen is sorted by drug name, changing a drug name to start with a non-Latin character causes the drug to display twice. When the issue occurs, editing either of the duplicate drug entries on the MDD properly changes the other entry. Re-sorting the MDD list or refreshing the MDD by switching to another page in the AT and then back causes the drug name to only display once, sorted after the names that start with Latin characters. If the drug is already in the SLS PCS formulary, then the drug name will display with the changed characters in the formulary and will be listed only once.
- ◆ **WARNING: The pre-import report does not identify duplicate drug IDs when importing via a CSV.** You are required to run a formulary report before generating an approved formulary and should review the Duplicate reports in the Safety worksheet to make sure all duplicates or potential duplicates in the MDD and formulary are identified.
- ◆ **In 2.6.x AT, user's cannot login with non-capitalized Active Directory groups.** In order to login, Active Directory Group names need to be capitalized when entered into the AT, otherwise Active Directory will not be usable.

### Known Uncommon Issues

This section details uncommon issues with 2.6.1 AT software that are unlikely to be experienced by most users. These issues are specific to the current release; previous issues can be found in previous Release Notes. None of the issues included in the document are considered to render the product unsafe or ineffective for clinical use.

- ◆ **WARNING: A concentration length of greater than 15 characters is converted to an incorrect value in the MDD when imported.** For example, a number such as 12345678901234567890 is converted to 92,233,720,368,547.77.
- ◆ **WARNING: If a drug record is modified in the MDD after it has been added to the formulary, the AT does not automatically update the audio file association.** If the drug name is changed in the MDD, the associated audio file will not change, resulting in a possible incorrect association. As a work around, you can either remove and then add the drug to the formulary to trigger a new audio file association, or manually select the correct audio file from the formulary edit dialog box.

- ◆ **WARNING: Combination drugs with long drug names do not display the second drug name. Combination drugs with very long drug names only display the first drug name as truncated.** There is no indication as to what the second drug name is. Users should look to modify the entry of the drug name (e.g. Calcium can be entered as Ca). Another option is to enter the drug name in the Tallman field of the AT which will reduce the font size.
- ◆ **Changing barcode parsing rule presets from US\_NDC to Custom does not warn you that custom rules will be removed.**
- ◆ **Zooming potentially causes a minor row offset in the AT.** This issue occurs only at 110% zoom in Chrome, and Edge (Chromium based). It does not occur in IE 11, Firefox, or older Edge. The AT actions column rows become slightly offset from the rest of the grid in all grids of the AT.
- ◆ **WARNING: Dilution units that are different from the parenteral container do not always provide a warning.** When the initial parenteral unit (e.g., mg) is the same as the first unit in the dilution (e.g., mg/mL) a warning is not provided. However, when the parenteral drug unit is different than the initial dilution units (e.g., parenteral container is g -> diluted container is mg/mL), the AT does provide a warning that they are different units.
- ◆ **Adding invalid feature keys displays an Error Contact Server (ECS).** When adding a feature key that has any of the following properties (e.g., empty field, single or double or special characters, all numbers or mix of numbers and characters) results in the AT displaying an ECS instead of Invalid Feature Key Error message.
- ◆ **SLS can lock-up when updating with a software app package (e.g., SLS-APP-2\_6\_0\_RC6-001.pkg).** When installing a multi-package on the SLS that contains a software APP package, the SLS can become stuck on the dialog indicating that the SLS was shutting down. If this occurs, power off and restart the SLS PCS, and reinstall the multi-package.
- ◆ **Second preparer's initials are not supported on a small label template pack (e.g., STD-S1 and ZH-S1).** Small labels do not support preparer's initials.

### Technical Support

If problems occur during software installation or operation, contact Codonics Technical Support at any time.

Phone: +1.440.243.1198  
 Email: support@codonics.com  
 Website: www.codonics.com

*Get it all with just one call*  
**800.444.1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.  
 Patents: www.codonics.com/ip/patents.

Copyright © 2010-2021 Codonics, Inc. Printed in the U.S.A. Part No. 901-245-026.02



17991 Englewood Drive Cleveland, OH 44130 USA +1 (440) 243-1198 +1 (440) 243-1334 Fax Email info@codonics.com www.codonics.com	Codonics Trading Co, Ltd. 317 Xianxia Rd. Building B Unit 1412 Changning Dist., Shanghai P.R. China, 200051 86-21-62787701 86-21-62787719 Fax	Codonics Limited KK AQUACITY 9F, 4-16-23, Shibaura Minato-ku, Tokyo, 108-0023 JAPAN 81-3-5730-2297 81-3-5730-2295 Fax
---	---	---