

# Codonics Horizon Operating Software

**Release Notes**

---

## **Version 2.1.1 Release Notes**

### **Summary**

Version 2.1.1 software provides support for a new type of Smart Card for all Codonics Horizon® Imagers.

The Smart Card is the component that provides licensing, configuration and feature information to the Horizon Imager.

The most recent domestic software version for Horizon Imagers was 2.0.0. The most recent international software version for Horizon Imagers was 2.1.0. Version 2.1.1 software unifies all Codonics Horizon Imagers to use one common software version along with the new type of Smart Card.

### **Products Affected**

This software upgrade and new Smart Card applies to all newly manufactured Horizon Imagers.

The Imager Software Version can be determined using the Control Panel by viewing the Imager Info screen from the Utilities menu.

### **Technical Support**

If problems occur during software installation, contact Codonics Technical Support between the hours of 8:30AM and 5:30PM EST (weekends and U.S. holidays excluded).

Phone: +1-440-243-1198

Email: [support@codonics.com](mailto:support@codonics.com)

Website: [www.codonics.com](http://www.codonics.com)

---

*Get it all with just one call*  
**1-800-444-1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2011 by Codonics, Inc. Printed in the U.S.A. Part No. 901-084-017.01



*We bring the future into focus*

17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)

Codonics Limited KK  
New Shibaura Bldg. F1  
1-3-11, Shibaura  
Minato-ku, Tokyo, 105-0023 JAPAN  
Phone: 81-3-5730-2297  
Fax: 81-3-5730-2295