

Studies show Codonics Safe Label System® is proven to increase patient safety

55% to ZERO


When clinicians prepared syringes manually in the OR, 55% contained a labeling error. When SLS was used in the OR, there were zero labeling errors and 100% compliance.

The operating room (OR) is a place where patients undergo complex and critical care procedures, under the vigilant eyes of talented doctors and nurses who are saving lives and curing disease. Despite vigilance and training, preventable drug administration errors continue to occur.

The Institute of Safe Medication Practices (ISMP) has identified many causes of medication errors, including incorrect dosing, accidental substitutions, and several other human errors, such as poor handwriting. Manual syringe labeling is often incomplete, incorrect or illegible, resulting in poor compliance with The Joint Commission (TJC) requirements and American Society of Anesthesiologists (ASA) guidelines, and creates a significant safety risk to patients.

Codonics Safe Label System (SLS) addresses the three most common medication errors made in the OR¹ and ensures medication safety and labeling compliance in the OR or anywhere medications are prepared. Every medication label prepared by SLS contains unique information for visual and 2D barcode-enabled BPOC (barcode point-of-care) confirmation of the prepared medication. Once a drug is prepared and labeled, the label allows for an additional visual and audible confirmation of the drug, concentration and beyond use date (BUD) pre-administration verification at SLS. The label also integrates with AIMS (Anesthesia Information Management Systems), other clinical information systems (CIS) and electronic health records (EHR).

Massachusetts General Hospital



Studies show SLS is proven to improve medication safety and labeling compliance

Massachusetts General Hospital (MGH) conducted a study in 2009 that looked at errors associated with clinician-prepared, manually labeled syringes. The study found that **55% of the syringes prepared manually by MGH clinicians contained a labeling error**, and the data was presented at the 2009 American Society of Anesthesiologists (ASA) Annual Meeting. In 2013, MGH anesthesiologists set out to assess the rates of perioperative medication errors and adverse drug events.² That study was published in the 2015 journal of Anesthesiology and reported that **zero syringes prepared using SLS contained errors.**³⁻⁴



SLS is an award-winning, FDA-cleared medical device, proven to:

- Help eliminate medication errors
- Improve patient safety
- Ensure consistency
- Provide 100% compliance every time a medication is prepared




In 2014, the University of Washington conducted a quality assurance audit whereby providers set out to determine the level of syringe labeling compliance among its anesthesia providers.⁵ The site randomly selected two cardiothoracic ORs and an OR satellite pharmacy over a five-week period for the pre-installation audit and over seven weeks for the post-installation audit, and both were conducted without the knowledge of the anesthesia providers. **After SLS was introduced, they found that 100% of the syringe labels prepared by SLS systems, whether by the anesthesia provider or satellite pharmacy, were compliant.** The study, published in the 2015 journal of Anesthesia & Analgesia, confirmed and concluded that **SLS significantly improves medication labeling compliance which can increase patient safety.**

The results of these studies showing the real-world impact of SLS establish that SLS is a “leapfrog safety innovation,” the benefit of which can be extended to every location where critical medications are being prepared by front line caregivers. SLS is a tool that is helping clinicians achieve “excellence every day.”

Today a standard of care in the world’s leading hospitals, SLS scans the FDA-mandated manufacturer’s vial/ampoule barcode at the point of care to create a medication label that complies with TJC requirements. SLS combines the FDA’s “source of truth” with your hospital’s rules for use and best practices in a drug database that is accessed on-demand every time a drug is scanned. Labels include the drug name, concentration, diluent and dilution, date and time of preparation and expiration, preparer’s initials and any warnings about the drug. The system also provides safety features including audio and visual feedback and clinical alerts for recalled and expired vials.

University of Washington

Post SLS installation compliance.⁵
Percent of prepared syringes in each group displayed by compliance and the number of label elements missing.



SLS has been used in more than 6 million drug preps, preventing over 45,000 medication errors, hundreds that would have resulted in death or serious injury.

“During a routine case, I needed Zofran. It wasn’t until I scanned the vial on Safe Label System that I realized I had Hydromorphone – a narcotic – which had been incorrectly placed in the anesthesia cart. Had it not been for Safe Label System’s visual and audible safety checks alerting me, my pediatrics patient may have never woke up.”

---Anesthesia Resident

1. Beverley A. Orser, Robert J.B. Chen and Doreen A. Yee. “Medication Errors in Anesthetic Practice: a Survey of 687 Practitioners.” Canadian Journal of Anesthesia 2001;48: 139-146

2. Nanji KC, Patel A, Shaikh S, et al. Evaluation of Perioperative Medication Errors and Adverse Drug Events. Anesthesiology 2015.

3. Vernest K, Nanji KC, Driscoll WD, et al. “Smart Labels”: Improving labeling compliance and patient safety in the operating room. American Society of Anesthesiologists Meeting 2009; Abstract A609.

4. Nanji KC, Vernest K, Sims N, et al. “Smart Labels”: Bar Code-Assisted Medication Labeling To Improve Efficiency and Patient Safety. American Society of Anesthesiologists Annual Meeting 2009. Accessed 2009, Abstract A612.

5. Jelacic, Srdjan MD; Bowdle, Andrew MD, et al. A System for Anesthesia Drug Administration Using Barcode Technology: The Codonics Safe Label System and Smart Anesthesia Manager™ Anesthesia & Analgesia 2015.

6. Based on Merry AF, Peck DJ. Anesthetists, errors in drug administration and the law. N Z Med J. 1995; 24:185-187

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patent Nos. 8,893,970; 8,639,525; 8,636,202; Patent(s) Pending. Copyright © 2010-2016 Codonics Inc. 03/2016



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com