


Safe Work Flow

- Improved medication safety with the Space Infusion System & Safe Label System™
- Helps eliminate drug and syringe swaps
- Improved accuracy to reduce typos
- Seamless integration using machine-readable barcode on label



2

All required information is included on syringe label
Syringe label barcode contains Drug Name and Concentration information



4

Scan label barcode and complete dose information
Helps eliminate syringe swaps



3

Space Infusion System

Insert labeled syringe in Perfusor Space



5

Drug is properly identified and Perfusor Space is ready to start
More accurate – less typos

Integration is work-in-progress

1

Safe Label System™

Scan unit dose container & automatically create a complete label

Helps eliminate drugs swaps