

# Codonics SLS 500i Operating Software

## Release Notes

### Version 1.8.2 SLS 500i Release Notes

#### Summary

Version 1.8.2 software is a software release for the Codonics Safe Label System™ SLS 500i™ Point of Care Station (PCS). Version 1.8.1 and 1.8.0 Administration Tool (AT) software can be used with the 1.8.2 PCS software. SLS PCS does not support software downgrades to previous software versions.

This document provides information about new features in 1.8.2 SLS PCS software. The SLS PCS User's Manual version 1.3.0 provides further information about SLS PCS.

**WARNING:** The use of label colors is intended only as an aid in the identification of drug groups and does not absolve the user from the duty of reading the label and correctly identifying the drug prior to use.

#### New Features

- ◆ **SLS supports a new Wi-Fi adapter.** A new Edimax Wi-Fi adapter provides support at both 2.4 and 5 GHz.

**NOTE:** Wi-Fi requires both an adapter and a feature key on SLS PCS. Contact your Codonics Sales Representative to obtain the proper adapter and feature key.

**NOTE:** Edimax Wi-Fi adapters can only be used with 1.8.2 or newer software. Contact Codonics Technical Support to upgrade existing SLS units using EnGenius and/or Asus Wi-Fi adapters to 1.8.2 or newer software.

#### Issues

- ◆ **Different configuration files may be required depending on the network security protocol used.** When using both the new Edimax Wi-Fi adapter and the existing EnGenius and/or Asus adapters, and using mixed mode WPA2/WPA Enterprise Network with certificates, contact Codonics Technical Support for specific instructions on setting up the configuration files.

#### Technical Support

If problems occur during software installation or operation, contact Codonics Technical Support at any time.

Phone: +1.440.243.1198

Email: support@codonics.com

Website: www.codonics.com

*Get it all with just one call*  
**1-800-444-1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2010-2014 Codonics, Inc. Printed in the U.S.A. Part No. 901-243-014.03.



17991 Englewood Drive  
Middleburg Heights, OH  
44130 USA  
+1.440.243.1198  
+1.440.243.1334 Fax  
Email info@codonics.com  
www.codonics.com

Codonics Trading Co, Ltd.  
317 Xianxia Rd. Building B  
Unit 1412  
Changning Dist., Shanghai  
P.R. China, 200051  
86-21-62787701  
86-21-62787719 Fax

Codonics Limited KK  
AQUACITY 9F,  
4-16-23, Shibaura  
Minato-ku, Tokyo,  
108-0023 JAPAN  
81-3-5730-2297  
81-3-5730-2295 Fax