



SAFETY


Scan

- Reads NDC barcode of FDA-approved drugs
- Visual and audible verification identifies and validates drug in hand
- Clinical alerts for recalled medications and expired syringes
- Prompts for dilutions* at preparation according to hospital rules

COMPLIANCE

Print

- Full color, easy-to-read labels on demand
- Complies with label requirements of TJC, APFS, USP and ASA guidelines
- Prints labels for IV lines and other OR needs



CONSISTENCY

Apply


- Labels include drug name, concentration, diluent/dilution*, date and time of preparation and expiration, preparer's initials and message*
- Easy-to-read and ready-to-apply to the prepared medication

AUTOMATIC DOCUMENTATION

Integrate

- Integrates with Epic and Cerner to record the exact NDC for improved charge capture, billing and 340B accuracy
- Less clicks for the provider reduces manual tasks
- BCMA in the OR

Exact NDC matching vial to syringe



Call Codonics today at **800.444.1198** or visit www.codonics.com for more information.

* When needed. All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: www.codonics.com/ip/patents Copyright © 2010-2020 Codonics Inc. 03/2020.